Stick®STEPPER / Stick®CARRIER

Cleaning of hand instruments

Stick®STEPPER**PROD 100035****Stick®CARRIER****PROD 100036**

ENGLISH

CLEANING OF Stick®STEPPER and Stick®CARRIER HAND INSTRUMENTS

StickSTEPPER and StickCARRIER instruments are specially designed to facilitate the use of Stick and everStick fibre products. The use of the instruments in different indications is described in the instructions for use of Stick and everStick fibre products.

DISINFECTION:

Recommendation for chemical solution:
It is recommended to use a pH 4,5–10,5 for silicone-handled instruments.

1. RINSING

Rinse the instrument with warm running water.

2. CLEANING ALTERNATIVES (A, B or C)

A. Cleaning by hand

- Remove visible impurities with warm water and a soft brush
- Use neutral detergent (pH ca. 7)
- Rinse with warm running water

B. Cleaning with an ultrasonic cleaner

- To protect instruments from wear, always have them in cassettes when cleaned ultrasonically
- After ultrasonic cleaning, rinse with warm running water

C. Washing machine

- Use neutral detergent (pH ca. 7)

3. STERILISATION ALTERNATIVES (A or B)

A. Autoclaving

- Autoclave clean, properly packed instruments

B. Dry-heat sterilisation

- Place dry, clean instruments in the dry-heat oven
- Do not load the oven too tightly
- The maximum temperature is +180°C (+355°F)

IMPORTANT:

- Have the machines serviced regularly, and follow the instructions for use; pay special attention to the correct dosage of chemical solution and exposure time
- Do not harm instruments by using hard brushes or engraved marks for identification, as these may crack the instrument surface

NOTE: Keep the instruments in the cassettes throughout the maintenance cycle.

WARRANTY: Stick Tech Ltd replaces defective products. Stick Tech does not assume liability for any damage or loss, whether direct or consequential, caused by improper use of the product or use contrary to the accompanying instructions. Before starting to use the product, the user is responsible for assessing the suitability of the

product for the intended purpose. The user bears all risk and liability for any consequences due to use contrary to the instructions or in inappropriate circumstances.

Manufactured for Stick Tech Ltd by LM Instruments.
For further information, please visit www.lminstruments.com.

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Fabriqué pour Stick Tech Ltd par LM Instruments.
Tillverkad för Stick Tech Ltd av LM Instruments.
Hergestellt von LM Instruments für Stick Tech Ltd.